NEW MODELS

SARTORIUS

Cubis® II

Milligram Balances



Standard Features

Modular design

Capacity: 150g-5,200g Readability: 0.0001g-0.01g

- User configurable and custom designed to meet changing applications
- Meets advanced Pharma compliance including user/password management for security, audit trail function logs and integrated alibi memory
- RS-232 and USB communication
- Easy cleaning by removing glass doors
- Underweigh hook included
- · Manual glass draft shield
- · Choice of display
- MCA Large, high-end, 7" color touchscreen TFT display in 16:9 format
- MCE State-of-the-art, TFT touchscreen operation

Options:

Bluetooth® wireless technology

Specifications

-p									
	3203S	2203S	2203P	1203S	5203S	5203P	623S	623P	323S
Part# List Price	195703	195704	195806	195705	195706	195807	195707	195808	195709
List Price	Consult	Consult	Consult	Consult	Consult	Consult	Consult	Consult	Consult
Part#	195711	195712	195809	195713	195714	195810	195715	195811	195716
Part# List Price	Consult	Consult	Consult	Consult	Consult	Consult	Consult	Consult	Consult
Capacity	3,200 g	2,200 g	1,010/2,200 g	1,200 g	5,200 g	5,200/2,400/1,200 g	620 g	150/300/620 g	320 g
Resolution	1 mg	1 mg	1/10 mg	1 mg	1 mg	5/2/1 mg	1 mg	1, 2, 5 mg	1 mg
Weighing Modes	Grams, Kilog	Grams, Kilograms, Carat, Pounds, Ounces, Troy Ounce, Tael Hong Kong, Tael Singapore, Tael Taiwan, Grain, Pennyweight, Milligram, Parts per Pound, Tael China, Mommes, Austrian Carat, Tola, Baht, Mesghal, Ton, Pound:Ounces, Newton, Micrograms							
Operating Modes (*)	Mass Unit Conversion, SQ Min Function, ISO-Cal Function, Density Determination (Buoyancy method only with MSE), Calculation, Averaging, Formulation, Percent Weighing, Totalizing, Counting, Customized Identification, Statistics, Time Control Function, Second Tare Memory, Over/Under Checkweighing								
Repeatability	1/6 mg	1 mg	1/6 mg	0.7 mg	1 mg	1 mg	0.7 mg	1/2/4 mg	0.7 mg
Linearity	5 mg	3 mg	5 mg	2 mg	5 mg	5 mg	2 mg	5 mg	2 mg
Stabilization Time	1 second	1 second	1 second	1 second	1 second	1 second	0.8 seconds	0.8 seconds	0.8 second
Operating Temperature		41° F to 104° F (5° C to 40° C)							
Weighing Chamber (DE)		(W x D x H) 8 x 8 x 7 in (from pan) (200 x 200 x 179 mm)							
Pan Size(W x D)		5.5 x 5.5 in (140 x 140 mm)							
Dimensions		(W x D x H) 16 x 9.8 x 11 in (410 x 249 x 281.4 mm) (with manual draft shield)							
Weight (Approx.)		17 lb							
Est. Ship Weight				27 lb				25 lb	
Calibration Weight	2,000 g	2,000 g	1,000 g	1,000 g	500 g	500 g	500 g	500 g	200 g
Power				120V	AC adapter (a	auto ranging)			
Warranty		Two-year Limited Warranty							
Return Policy		Sartorius Standard Policy - restock fee may apply							
*MCF Large TET Touchscreen a	nd agev_to_uea ma	nu evetom							

^{*}MCE Large TFT Touchscreen and easy-to-use menu system

^{**}MCA TFT color touchscreen, Ethernet and SD card reader

NEW MODELS SARTORIUS Cubis® II Milligram Balances

Part Number/Price

Part#	Model#	Description	Price
195703*	MCE3203S-2S00-E	3,200 g x 1 mg	Consult
195704*	MCE2203S-2S00-E	2,200 g x 1 mg	Consult
195806*	MCE2203P-2S00-E	1,010 g/2,200 g x 1 mg/10 mg	Consult
195705*	MCE1203S-2S00-E	1,200 g x 1 mg	Consult
195706	MCE5203S-2S00-E	5,200 g x 0.001 g	Consult
195807	MCE5203P-2S00-E	5,200 g/2,400 g/1,200 g x 0.005 g/0.002 g/0.001 g	Consult
195707*	MCE623S-2S00-E	620 g x 1 mg	Consult
195808*	MCE623P-2S00-E	150 g/300 g/620 g x 1 mg/2 mg/5 mg	Consult
195709*	MCE323S-2S00-E	320 g x 1 mg	Consult
195711**	MCA3203S-2S00-E	3,200 g x 1 mg	Consult
195712 **	MCA2203S-2S00-E	2,200 g x 1 mg	Consult
195809 **	MCA2203P-2S00-E	1,010 g/2,200 g x 1 mg/10 mg	Consult
195713 **	MCA1203S-2S00-E	1,200 g x 1 mg	Consult
195714	MCA5203S-2S00-E	5,200 g x 0.001 g	Consult
195810	MCA5203P-2S00-E	5,200 g/2,400 g/1,200 g x 0.005 g/0.002 g/0.001 g	Consult
195715 **	MCA623S-2S00-E	620 g x 1 mg	Consult
195811 **	MCA623P-2S00-E	150 g/300 g/620 g x 1 mg/2 mg/5 mg	Consult
195716 **	MCA323S-2S00-E	320 g x 1 mg	Consult

^{*}MCE Large TFT Touchscreen and easy-to-use menu system

Options/Accessories

Part#	Model#	Description	Price
Consult	TR	Certificate of compliance and full factory test report (factory installed)	Consult
Consult	IR	Additional RS-232, 25-pin connection (factory installed)	Consult
156906	6971987	Power supply, 100-240V	Consult
156908	6900900	Cord, power 220V EU	Consult
Consult	IB	Bluetooth® wireless interface (factory installed)	Consult
Consult	IP	Additional RS-232, nine-pin connection with PS2 port (factory installed)	Consult
Consult	DA, DR, DI or DU	Other draft shields available	Consult

^{**}TFT color touchscreen, Ethernet and SD card reader